GB Instructions for cleaning and sterilization of stainless steel instrumentation

Operation conditions
Temperature 10° C to 30° C. **Humidity** 30% to 85% R.H.
Storage/Transportation Conditions
 The packing of the instruments must be able to withstand the transportation conditions and protect the instruments from damage. Store in a cool dry place, away from direct sunlight.
Temperature 05° C to +50° C.

General recommendations
 Reusable instruments are provided **NON-STERILE** and must be cleaned and sterilized according to these instructions prior to use.

Manual Cleaning
Equipment: Latex Gloves, proteolytic enzyme, non-metallic brush, running water, purified water, container for solution, container for rinsing.
Step 1: Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
Step 2: Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulas should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
Step 3: Soak instruments for a minimum of 10 minutes. While soaking scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces and areas with moving components or springs. Lumens, blind holes and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round brush into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.
Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.
Step 4: Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
Step 5: Dry instruments with a clean, absorbent non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas.

Combination Manual/Automated Cleaning Steps
 Steps 1 to 4: same as above
Step 5: Place instruments in a suitable validated washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure; e.g. open all instruments, place concave instruments on their side or upside up, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfector is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer's instructions. According to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash • Cold Softened Tap Water • 2 minutes
2	Enzyme Spray & Soak • Hot Softened Tap Water • 1 minute
3	Rinse • Cold Softened Tap Water
4	Detergent Wash • Hot Tap Water (64-66 °C) • 2 min.
5	Rinse • Hot Purified Water (64-66 °C) • 1 minute
6	Hot Air Dry (116 °C) • 7 to 30 minutes

Note: The washer/disinfector manufacturer's instructions should be followed.

Sterilization
 Most heat/steam sterilization is recommended for reusable stainless-steel instruments. Use of an approved chemical integrator or chemical emulater within each sterilization load is recommended. Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. The manufacturer's recommendations for installation, validation, and maintenance should be followed. Validated exposure times and temperatures to achieve a 10-6 sterility assurance level (SAL) are:

Cycle type	Minimum temperature	Minimum exposure time	Pressure
Pre-vacuum / Vacuum Pulse	134 °C	3 minutes	28 p.s.i.

The recommended drying time for single wrapped instruments is 20 minutes unless otherwise noted in device specific instructions. A 30 minutes minimum cooling time is recommended after drying but longer times may be necessary because of load configuration, ambient temperature and humidity, device design and packaging used.
Notes: Disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern about TSE/CJD contamination are: 134 °C/273 °F for 18 minutes. Our stainless-steel instruments are compatible with these parameters.

Limitations on Reprocessing
 Repeated processing according to these instructions has minimal effect upon metal reusable instruments and accessories unless otherwise noted. End of life for stainless steel or other metal instruments is generally determined by wear and damage incurred during the intended use. Reusable instruments comprised of polymers or incorporating polymer components can be sterilized using steam however they are not as durable as their metal counterparts. If polymer surfaces show signs of excessive surface damage (e.g. crazing, cracks or delamination), distortion or are visibly warped they should be replaced. Non-foaming, neutral pH enzymatic and cleaning agents are recommended for processing reusable instruments and accessories. Alkaline agents with a pH of 12 or less may be used to clean stainless steel and polymer instruments in countries where required by law or local ordinance or where prior diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt-Jakob Disease (CJD) are a concern. It is critical that alkaline cleaning agents are completely and thoroughly neutralized.

 Medical Device

 Manufacturer

 EU Representative

 Catalogue #

 Lot Code

 Date of Manufacture

 CE Marked

Notice: In the event of a serious incident occurring in connection with this instrument, the user or patient of this device shall report directly to the manufacturer and the competent authority where the user of patient is established.

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